

## DispensMate Bottle-Top Dispenser



### 제품특징

- Excellent chemical resistance, components are made of PTFE, FEP, BSG, PP
- Fully autoclavable at 121°C
- Four ranges of bottle-top dispenser covering a volume range from 0.5mL to 50mL
- Easy for cleaning and maintenance
- The optional flexible discharge tube with safety handle permits fast and precise dispensing
- Vapor pressure Max. 500mbar, viscosity max. 500mm<sup>2</sup>/s, temperature max. 40°C, density max. 2.2g/cm<sup>3</sup>
- Dispensmate is supplied with S40, GL32, GL38, GL25, GL28 sized adaptors

### Accessories

- 17100553 GL45, Adapter, 45mm
- 17000069 S40, Adapter, 45/40mm
- 17400017 GL32, Adapter, 45/32mm
- 17400018 GL38, Adapter, 45/38mm
- 17400019 GL25, Adapter, 32/25mm
- 17400020 GL28, Adapter, 32/28mm
- 17400037 Reagent Bottle (Brown, 1L, Options)
- 17400021 Filling tube(m)
- 17100073 Discharge tube(m)

### 제품사양

Cat. No.	Volume Range (mL)	Graduation(mL)	A <sub>s</sub> ±		CV≤	
			%	μL	%	μL
731100017777	0.5~5	0.1	0.5	25	0.1	5
731100027777	1.0~10	0.2	0.5	50	0.1	10
731100037777	2.5~25	0.5	0.5	125	0.1	25
731100047777	5.0~50	1.0	0.5	250	0.1	50