


## StepMate Stepper



### 제품특징

- Light and ergonomic solution provide easy usability
- Single-handed operation
- Volume range from 1 $\mu$ L to 5mL
- Up to 48 dispensing steps
- Maintenance free
- Equipped with durable tip
- insertion lever
- Fits seven sizes of disposable polypropylene syringes from 0.5mL to 50mL

### How to select dispensing volume

- Find the required dispensing volume in the table
- Use the adjusting wheel to set the required volume
- Select and insert suitable tips
- Follow the maximum dispensing steps

### 제품사양

Cat. No.	Syringe Volume(mL)	Sample volume					A %	P %	Description
17900032	0.5	10	20	30	40	50	$\pm 0.8$	$\leq 0.7$	0,5 ml, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900033	1.25	25	50	75	100	125	$\pm 0.8$	$\leq 0.5$	1,25 ml, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900034	2.5	50	100	150	200	250	$\pm 0.8$	$\leq 0.4$	2,5 ml, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900035	5	100	200	300	400	500	$\pm 0.4$	$\leq 0.3$	5ml, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900036	12.5	250	500	750	1000	1250	$\pm 0.3$	$\leq 0.3$	12,5 ml, piston PE-HD, cylinder PP, 1 packing, 5 pcs/packing
17900037	25	500	1000	1500	2000	2500	$\pm 0.2$	$\leq 0.3$	25ml, piston PE-HD, cylinder PP, 1 packing, 2pcs/packing with 1 pc adaptor
17900038	50	1000	2000	3000	4000	5000	$\pm 0.2$	$\leq 0.2$	50ml, piston PE-HD, cylinder PP, 1 packing, 2pcs/packing with 1 pc adaptor