

## FlatSpin Ultra-flat Compact Magnetic Stirrer

With excellent stirring functionality, the compact shape of FlatSpin takes less bench space and makes storage easier

### 제품특징

- A wide range of speed from 15 to 1500rpm
- Automatic reverse rotation for excellent mixing
- Protection class (IP 65)
- PET top provides excellent chemical resistance
- Ultra silence and vibration-free



### 제품사양

Model	FlatSpin
Dimension of work plate	φ90mm
Work plate material	PET
Motor type	Motorless
Motor output power	3W
Power	5W
Voltage	100-120/200-240V, 50/60Hz
Stirring positions	1
Max. stirring quantity [H2O]	0.8L
Max. magnetic bar [length]	40mm
Speed range	15-1500rpm
Speed setting	button
Protection class	IP65
Dimension [W x D x H]	175×120×15mm
Weight	0.3kg
Permissible ambient temperature and humidity	5-40°C, 80%