

KSS-750F / KSS-1500F / KSS-2500F Table Sonicators 초음파 추출.파쇄.균질.분산.유화기



Big Volume Power Adjustable  
Nano Sonics Ultrasonic Homogenizer

재료,소재 분산을 위해 특별히 설계된 초음파 처리 장비로 강력한 출력과 우수한 분산 성능을 보여줍니다

제품특징

- With new software central microcomputer centralized control
- It adopts Backlit and LCD large screen display
- Ultrasonic power step with 1% continuous adjustable.
- With pulse and continuous operation, and with test function
- Ultrasonic time, gap time, and total time can be accurate to 0.1 seconds.

제품사양

Model	KSS-750F	KSS-1500F	KSS-2500F
Frequency	19.5-20.5 KHz	19.5-20.5 KHz	19.5-20.5 KHz
Display	LCD	LCD	LCD
Adjustable power range	750W (20%-99%)	1500W (20%-99%)	2400W (20%-99%)
Accompany Tip	10mm	20mm	38mm
Optional Tip	6,15 mm	15, 25 mm	None
Processing Volume	5-700 ml	50-1200 ml	500-5000 ml
Duty Ratio	0.1-99.9%		
Temperature Alarm	Yes(Prevent sample over heating)		
Alarm	time, temperature, overload		
The size of Power Supply Chassis	400x280x220(mm)	400x280x220(mm)	400x280x220(mm)
Standard configuration	Generator host 1pc	Generator host 1pc	Generator host 1pc
	Dispersing head 1pc	Dispersing head 1pc	Dispersing head 1pc
	Sound Abating Chamber 1pc	Sound Abating Chamber 1pc	Sound Abating Chamber 1pc