

MS-C-S4 4-position Magnetic Stirrer for Cell Culture

MS-C series provide remarkable performance in constant speed stirring, cyclic stirring and dual speed stirring with various working modes for cell culture.

제품특징

- Programmable with multi mode selection
- Digital LED displays the actual speed
- Available in both single position and 4 position
- Low speed prevents cell damage



제품사양

Model	MS-C-S1
Dimension of work plate	φ105mm
Work plate material	Cast Aluminum
Motor type	Stepper Motor
Motor output power	15W
Voltage	100~120/200~240V, 50/60Hz
Stirring positions	4 (each point distance 190mm)
Max. stirring quantity [H ₂ O]	3L
Max. magnetic bar [length]	80mm
Speed range	5-200rpm
Speed display	LED
Timer function	0-49min59s
Data connector	RS232
Protection class	IP21
Dimension [W × D × H]	395 × 450 × 90mm
Weight	5.5kg
Permissible ambient temperature and humidity	5-40°C, 80%RH