

MX-T6-S Tube Roller

DLAB tube rollers provide a gentle and highly efficient rolling motions, ideal for mixing blood samples, viscous substances and liquid-solid suspensions. Used in a variety of applications, such as the prevention of blood coagulation, immune precipitation, etc.

제품특징

- Speed range of 0~70rpm
- Six roller design for gentle rolling and easy to clean



제품사양

Model	MX-T6-S
Motor type	DC motor
Rocking motion	Rocking and rolling
Amplitude	24mm
Max. load capacity	4kg
Number of rollers	6
Roller size [length]	280mm
Speed range	0-70rpm
Speed display	Scale
Operation mode	Continuous
Voltage	100-240V, 50/60Hz
Power	25W
Dimensions [D × W × H]	450 × 260 × 120mm
Weight	4.5kg
Permissible ambient temperature and humidity	5~40°C, 80%RH
Protection class	IP21