

SK-L180-Pro LCD Digital Linear Shaker

DLAB orbital and linear shakers are ideal for almost any vessel from microcentrifuge tubes through petri dishes and microtitre plates to conical flasks.

제품특징

- Max. load capacity is 7.5kg(L330) and 2.5kg(L180) with platform
- Double LCDs for independent display of time and speed
- Wide range of platforms for use with a variety of vessels
- Maintenance- free brushless DC motor
- Over speed detection and protection



제품사양

Model	SK-L180-PRO
Orbital diameter	4mm
Max. load capacity (with platform)	2.5kg
Motor type	Brushless DC motor
Motor input power	28W
Motor output power	15W
Speed range	100~350rpm
Speed display	LCD
Timer display	LCD
Time settings range	1min~19h59min
Voltage	100~240V,50/60Hz
Power	30W
Weight	8.1kg
Permissible ambient temperature and humidity	5~40°C,80%RH
Dimension[W×D×H]	300×340×100mm
Protection class	IP21
Data connector	RS232