

## SK-O180-S LED Digital Orbital Shaker

DLAB orbital and linear shakers are ideal for almost any vessel from microcentrifuge tubes through petri dishes and microtitre plates to conical flasks.

### 제품특징

- Max.load capacity is 3kg with platform
- Double LEDs for independent display of time and speed
- Speed range of 40-200rpm
- Wide range of platforms for use with a variety of vessels



### 제품사양

Model	SK-O180-S
Orbital diameter	20mm
Max. load capacity (with platform)	3kg
Motor type	DC motor
Motor input power	16W
Motor output power	10W
Speed range	40~200rpm
Speed display	LED
Timer display	LED
Time settings range	1min~19h59min
Voltage	100~240V,50/60Hz
Power	20W
Weight	3.1kg
Permissible ambient temperature and humidity	5~40°C,80%
Dimension[W×D×H]	270×330×110mm
Protection class	IP21
Data connector	-