

## SK-O330-Pro LCD Digital Orbital Shaker

DLAB orbital and linear shakers are ideal for almost any vessel from microcentrifuge tubes through petri dishes and microtitre plates to conical flasks.

### 제품특징

- Max. load capacity is 7.5kg(O330) and 2.5kg(O180) with platform
- Double LCDs for independent display of time and speed
- Used with a wide variety of accessories
- Maintenance- free brushless DC motor
- Over speed detection and protection
- Remote function provides PC control and data transmission



### 제품사양

Model	SK-O330-PRO
Orbital diameter	10mm
Max. load capacity (with platform)	7.5kg
Motor type	Brushless DC motor
Motor input power	28W
Motor output power	15W
Speed range	100~500rpm
Speed display	LCD
Timer display	LCD
Time settings range	1min~19h59min
Voltage	100~240V,50/60Hz
Power	30W
Weight	13.5kg
Permissible ambient temperature and humidity	5~40°C,80%
Dimension[W×D×H]	370×420×100mm
Protection class	IP21
Data connector	RS232