

SK-R1807-S LED Digital Rocking shaker

DLAB rocking and 3D shakers provide perfect solution in electrophoresis gel dyeing/dicoloration, sample washing, molecular hybridization, immunoprecipitation, Western blot and cell culture.

제품특징

- Rocking motion angle 7°
- Max. load capacity is 3kg with platform
- Speed range of 10-80rpm
- Double LEDs for independent display of time and speed
- Wide range of platforms for use with a variety of vessels



제품사양

Model	SK-R1807-S
Tilt angle	7°
Max. load capacity (with platform)	3kg
Motor type	DC motor
Motor input power	16W
Motor output power	10W
Speed range	10-80rpm
Speed display	LED
Timer display	LED
Time settings range	1min-19h59min
Voltage	100-240V, 50/60Hz
Power	20W
Weight	3.2kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Dimensions [W × D × H]	270 × 330 × 130mm
Protection class	IP21