

## T25 digital / DI25 digital 파쇄.균질.분산.유화기



### 제품특징

High-performance dispersing instrument for volumes from 1 to 2,000ml (H<sub>2</sub>O).  
 The spectrum of applications ranges from homogenizing waste water sample to the use in laboratory reactors, to dispersion tasks under vacuum / pressure and sample preparation in medical diagnostics.

- Three types of shaft bearings
- Standard version with digital display
- Rotor-Stator configurations have thirty years of proven, guaranteed comparability of test results.

### 제품사양

Model	T25 digital	DI25 digital
Ident No.	3565000	3717000
Motor capacity	500W	500W
Speed range	3,400~24,000 rpm	3,400~24,000 rpm
Speed display	Digital	Digital
Speed adjustment	Stepless	Stepless
Volume range (H <sub>2</sub> O)	1 ~ 2,000 ml	1 ~ 2,000 ml
Max. viscosity	5,000 mPas	5,000 mPas
Noise without dispersing element	73 dB	73 dB
Overload protection	Yes	Yes
Smooth run-on/start	Yes	Yes
Dimension [W x D x H]	65 x 80 x 240mm	65 x 80 x 240mm
Weight	1.6 kg	1.6 kg