

## D1008 series Palm Micro Centrifuge

### 제품특징


- Ideal for microfiltration and quick spin-downs
- Easy-to-use with dual start/stop function
- 7000rpm fixed speed(D1008)
- 5000rpm fixed speed(D1008E)
- Compatible with 2mLx8 and PCR8x2 rotors
- Quiet running, noise level  $\leq 45$ dB
- Rotor change without tools
- Accelerates and brakes in seconds
- Available in two colors i.e., green, blue



### 제품사양

Model	D1008	D1008E
Max. Speed	7000rpm	5000rpm
Max. RCF	2,680 xg	1360 xg
Rotor Capacity	0.2/0.5/1.5/2.0mL x 8; 0.2mL x 16 PCR strips or 0.2mL x 2 PCR 8 strips	0.2/0.5/1.5/2.0mL x 8; 0.2mL x 16 PCR strips or 0.2mL x 2 PCR 8 strips
Run Time	Continuous operation	Continuous operation
Driving Motor	DC motor	DC motor
Power	AC110-240V/50Hz/60Hz 20W	AC110-240V/50Hz/60Hz 20W
Noise level(dB)	$\leq 45$ dB	$\leq 45$ dB
Dimensions [W x D x H]	160x170x122mm	160x170x122mm
Weight	0.5kg	0.5kg

## Available rotors



Model	A08-2	A02-PCR8
Cat. No.	19400015	19400017
Max. Speed	7,000rpm	10000/7000/5000rpm
Max. RCF	F:5460/2680/1360×g	5460/2680/1360×g
Rotor Capacity	2mL/1.5mL×8 0.2mL×8,0.5mL×8	0.2mL×16 PCR strips, 0.2 mL×2 PCR 8 strips
Rotor material	High strength plastic	High strength plastic

## Available adapters

### ❖ SA02P2

- 0.2mL rotor adapter
- Used with A08-2 rotor
- 8pcs/pk

### ❖ SA05P2

- 0.5mL rotor adapter
- Used with A08-2 rotor
- 8pcs/pk