

## MS7-S 7" Square Magnetic Stirrer

### 제품특징

- Speed range of 0-1500rpm
- Max. stirring quantity of H2O at 10L
- Ceramic work plate provides excellent chemical resistant performance



### 제품사양

Model	MS7-S
Dimension of work plate	184×184mm ( 7inch )
Work plate material	Glass ceramic
Motor type	Shaded pole motor
Motor input power	15W
Motor output power	1.5W
Power	30W
Voltage	100-120/200-240V, 50/60Hz
Stirring positions	1
Max. stirring quantity [H2O]	10L
Max. magnetic bar [length]	80mm
Speed range	0-1500rpm
Speed display	Scale
Protection class	IP2X
Dimension [W×D×H]	215×360×112mm
Weight	3.8kg
Permissible ambient temperature and humidity	5-40°C, 80%