

## MS-M-S10 10-Position Magnetic Stirrer

### 제품특징

- Maintenance-free brushless DC motor
- Speed range of 0-1100rpm
- Stainless steel work plate covered with silicone cushion, provides excellent performance of skid resistance



### 제품사양

Model	MS-M-S10
Dimension of work plate[W x D]	180×450mm
Work plate material	Stainless steel with silicone coating
Motor type	Brushless DC motor
Motor input power	12W
Motor output power	4W
Power	20W
Voltage	100~120/200~240V,50/60Hz
Stirring positions	10
Max. stirring quantity [H <sub>2</sub> O]	0.4LX10
Max. magnetic bar [length]	40mm
Speed range	0~1100rpm
Speed display	Scale
Protection class	IP21
Dimension [W x D x H]	182×622×65mm
Weight	3.2kg
Permissible ambient temperature and humidity	5-40°C, 80%RH