

MS-S Magnetic Stirrer

제품특징

- Maintenance free brushless DC motor
- Speed range of 0-1500 rpm
- Max. stirring quantity of H₂O at 20L
- Stainless steel work plate cover with ceramic material provides good chemical resistance performance



제품사양

Model	MS-S
Dimension of work plate	φ135(5 inch)mm
Work plate material	Stainless steel cover with ceramic
Motor type	Brushless DC motor
Motor input power	18W
Motor output power	10W
Power	30W
Voltage	100-120/200-240V, 50/60Hz
Stirring positions	1
Max. stirring quantity [H ₂ O]	20L
Max. magnetic bar [length]	80mm
Speed range	0-1500rpm
Speed display	Scale
Protection class	IP2X
Dimension [W×D×H]	160×280×85mm
Weight	2.8kg
Permissible ambient temperature and humidity	5-40°C, 80%