

## SK-D1807-E Economical 3D Shaker

DLAB rocking and 3D shakers provide perfect solution in electrophoresis gel dyeing/dicoloration, sample washing, molecular hybridization, immunoprecipitation, Western blot and cell culture.

### 제품특징

- Three-dimensional rocking motion with the angle of 7°
- Max. load capacity is 3kg with platform
- Speed range of 0-80rpm
- Wide range of platforms for use with a variety of vessels



### 제품사양

Model	SK-D1807-E
Tilt angle	7°
Max. load capacity(with platform)	3kg
Motor type	DC motor
Motor input power	16W
Motor output power	10W
Speed range	0-80rpm
Speed display	Scale
Timer display	100-240V,50/60Hz
Time settings range	20W
Voltage	2.7kg
Power	5-40°C, 80%RH
Weight	270×330×130mm
Permissible ambient temperature and humidity	IP21
Dimensions[W × D × H]	270×330×130mm
Protection class	IP21